

# FillerCollector™



Gentle. Effortless. Efficient.

FillerCollector™ is a better way to harvest fat.

## Ordering Information

| Description                     |                         | Units   | REF      |
|---------------------------------|-------------------------|---|----------|
| FillerCollector™                |                         | Non-sterile, autoclavable Complete Set, including 2 Disposable Sets (sterile silicone washer and plug), 2 Extraction Cannulas | 660000   |
| Disposable Set FillerCollector™ | Sterile, for single use | 10 Sets   | 660023   |
|                                 |                         | 5 Sets  | 660023-5 |
| Extraction Cannula 300mm        | Sterile, for single use | 10 Sets   | 665010   |
|                                 |                         | 5 Sets  | 665010-5 |

Brochure not for distribution in the USA. FillerCollector™ is a trademark owned by human med AG, Germany. Patents pending.



©human med AG, 09/2010 REF 9001008

human med AG  
Wilhelm-Hennemann-Str. 9  
19061 Schwerin, Germany  
Tel: +49 (0) 385 / 3 95 70 -0  
Fax: +49 (0) 385 / 3 95 70 -10  
info@humanmed.com

# The FillerCollector™ is a better way to harvest fat.

In the past, fat harvesting was often a difficult, time-consuming process with unpredictable results.

But now, the all new FillerCollector™ simplifies the collection of injectable fat for autologous fat transfer. The FillerCollector™ is the most advanced fat harvesting system available for small to medium amounts of fat. The filterless, reusable collection device for fat harvesting, the FillerCollector™ is fast, easy and gentle on fat.

## Gentle

- Gentle on fat to help improve the chances of long-term fat survival\*
- No centrifuging or other processing needed - saves time and protects fat cells from harm
- Pressure equalization inside container helps protect fat cells
- Gentle, passive blending of fat within container - improves fat consistency from multiple harvest sites

## Effortless

- Simplified collection process
- Unique design continuously separates fat from rinsing fluid - no additional washing needed
- Easy fat extraction from container
- Filterless - no clogging

## Efficient

- Harvests high-quality fat\*\*, which is ready for immediate transfer
- Pre-filter traps fibrous strains for more homogeneous fat
- Sterile, closed loop system
- Robust and durable components
- Autoclavable
- Harvests 100 ml in just 10 minutes

\*P. Smith et al: Autologous Human Fat Grafting: Effect of Harvesting and Preparation Techniques in Adipocyte Graft Survival. Plastic and reconstructive surgery, 117: 1836-44 (2006)

\*\* Gordon H. Sasaski, M.D.: Fat Viability Study, Aesthetic Guide Nov/Dec 2009



For gentle and effortless fat harvesting