

Material needed for surgery



Liposuction	Fat transfer
<p><b>(D)</b> WAL Applicator (Ref.: 500001)</p> <p><b>(R)</b> Ø 2,5 mm Infiltration cannula (Ref.: 500070) or Ø 3,5 mm Infiltration cannula (Ref.: 500028)</p> <p><b>(R)</b> Ø 3,8 mm Rapid (Ref.: 501092), Ø 4,2 mm (Ref.: 500149) irrigation / aspiration WAL cannula</p>	<p><b>(D)</b> WAL Applicator (Ref.: 500001)</p> <p><b>(R)</b> Ø 3,5 mm STS Infiltration cannula (Ref.: 503003)</p> <p><b>(R)</b> Ø 3,8 mm STS Rapid irrigation / aspiration WAL cannula (Ref.: 503002)</p>
	<p>LipoCollector (large amount) (Ref.: 670000)</p> <ul style="list-style-type: none"> <li><b>(R)</b> Collection container, lid, basket, suction pipe, separation pipe, foot ring, water straw must be sterilized</li> <li><b>(D)</b> Disposable set (includes silicone sealing ring, sealing plug) (Ref.: 670200-5)</li> <li><b>(D)</b> Connection tubing (connection between collector and suction bag) (Ref.: 1306.3000.100)</li> <li><b>(D)</b> Extraction cannula (Ref.: 675010-5)</li> </ul> <p>FillerCollector (small amount) (Ref.: 660000)</p> <ul style="list-style-type: none"> <li><b>(R)</b> Collection container, lid, basket, suction pipe, separation pipe, foot ring, water straw must be sterilized</li> <li><b>(D)</b> Disposable set (includes silicone sealing ring, sealing plug) (Ref.: 660023-5)</li> <li><b>(D)</b> Connection tubing (connection between collector and suction bag) (Ref.: 1306.3000.100)</li> <li><b>(D)</b> Extraction cannula (Ref.: 665010-5)</li> </ul>
	<p><b>(D)</b> Luer-to-Luer connector (fem/fem) (Ref.: 5206634) to transfer the extracted fat from 50 cc syringes into smaller syringes</p>
	<p>Refill</p> <p>Large amount</p> <p><b>(D)</b> BEAULI® infiltration cannula (Ref.: 655020-5)</p> <p>Small amount</p> <p><b>(D)</b> BEAULI® infiltration cannula (Ref.: 655030-5)</p>
<p>Other</p> <ul style="list-style-type: none"> <li><b>(D)</b> 3000 ml suction bags (Ref.: 57187)</li> <li>If not connected or damaged Hydrophobic Filter-Set, for one side (Ref.: 101947)</li> </ul>	<p>Other</p> <ul style="list-style-type: none"> <li><b>(D)</b> 3000 ml suction bags (Ref.: 57187)</li> <li>If not connected or damaged Hydrophobic Filter-Set, for one side (Ref.: 101947)</li> <li><b>(R)</b> Syringe rack (Ref.: 610300) or other syringe holder</li> </ul>

**(D)** = Disposable  
**(R)** = Reusable