

## Material needed for surgery



Liposuction	Fat transfer
Small amount, e.g. chin (D) Biofill Applicator (Ref.: 580001) (R) Biofill cannula STS, Irrigation / Aspiration Cannula, Ø 3,5 mm (Ref.: 580086)	Small amount, e.g. scar correction (D) Biofill Applicator (Ref.: 580001) (R) Biofill cannula STS, Irrigation / Aspiration Cannula, Ø 3,5 mm (Ref.: 580086)
	FillerCollector (small amount) (Ref.: 660000) <ul style="list-style-type: none"> <li>(R) Collection container, lid, basket, suction pipe, separation pipe, foot ring, water straw must be sterilized</li> <li>(D) Disposable set (includes silicone sealing ring, sealing plug) (Ref.: 660023-5)</li> <li>(D) Connection tubing (connection between collector and suction bag) (Ref.: 1306.3000.100)</li> <li>(D) Extraction cannula (Ref.: 665010-5)</li> </ul>
	(D) Luer-to-Luer connector (fem/fem) (Ref.: 5206634) to transfer the extracted fat from big syringes into smaller syringes
	Refill (D) Small amount BEAULI® infiltration cannula (Ref.: 655030-5)
Other <ul style="list-style-type: none"> <li>(D) 1000 ml suction bags (Ref.: 57157)</li> <li>If not connected or damaged Hydrophobic Filter-Set, for one side (Ref.: 101947)</li> </ul>	Other <ul style="list-style-type: none"> <li>(D) 1000 ml suction bags (Ref.: 57157)</li> <li>If not connected or damaged Hydrophobic Filter-Set, for one side (Ref.: 101947)</li> <li>(R) Syringe rack (Ref.: 610300) or other syringe holder</li> </ul>

(D) = Disposable

(R) = Reusable