

body-jet® evo

Liposuction	Fat transfer
<p>Large amount</p> <p>(D) WAL Applicator, 370 cm (Ref.: 570001)</p> <p>(R) Ø 2,5 mm Infiltration cannula (Ref.: 500070) or Ø 3,5 mm Infiltration cannula (Ref.: 500028)</p> <p>(R) Ø 3,8 mm Rapid (Ref.: 501092), Ø 4,2 mm (Ref.: 500149) irrigation / aspiration WAL cannula</p>	<p>Large amount</p> <p>(D) WAL Applicator, 370 cm (Ref.: 570001)</p> <p>(R) Ø 3,5 mm STS Infiltration cannula (Ref.: 503003)</p> <p>(R) Ø 3,8 mm STS Rapid irrigation / aspiration WAL cannula (Ref.: 503002)</p>
<p>Small amount, e.g. chin</p> <p>(D) Biofill Applicator (Ref.: 580001)</p> <p>(R) Biofill cannula STS, Irrigation / Aspiration Cannula, Ø 3,5 mm (Ref.: 580086)</p>	<p>Small amount, e.g. scar correction</p> <p>(D) Biofill Applicator (Ref.: 580001)</p> <p>(R) Biofill cannula STS, Irrigation / Aspiration Cannula, Ø 3,5 mm (Ref.: 580086)</p>
	<p>LipoCollector (large amount) (Ref.: 670000)</p> <ul style="list-style-type: none"> (R) Collection container, lid, basket, suction pipe, separation pipe, foot ring, water straw must be sterilized (D) Disposable set (includes silicone sealing ring, sealing plug) (Ref.: 670200-5) (D) Connection tubing (connection between collector and suction bag) (Ref.: 1306.3000.100) (D) Extraction cannula (Ref.: 675010-5) <p>FillerCollector (small amount) (Ref.: 660000)</p> <ul style="list-style-type: none"> (R) Collection container, lid, basket, suction pipe, separation pipe, foot ring, water straw must be sterilized (D) Disposable set (includes silicone sealing ring, sealing plug) (Ref.: 660023-5) (D) Connection tubing (connection between collector and suction bag) (Ref.: 1306.3000.100) (D) Extraction cannula (Ref.: 665010-5)
	<p>(D) Luer-to-Luer connector (fem/fem) (Ref.: 5206634) to transfer the extracted fat from 50 cc syringes into smaller syringes</p>
	<p>Refill</p> <p>Large amount</p> <p>(D) BEAULI® infiltration cannula (Ref.: 655020-5)</p> <p>Small amount</p> <p>(D) BEAULI® infiltration cannula (Ref.: 655030-5)</p>
<p>Other</p> <ul style="list-style-type: none"> (D) 3000 ml suction bags (Ref.: 57187) If not connected or damaged Hydrophobic Filter-Set, for one side (Ref.: 103532) 	<p>Other</p> <ul style="list-style-type: none"> (D) 3000 ml suction bags (Ref.: 57187) If not connected or damaged Hydrophobic Filter-Set, for one side (Ref.: 103532) (R) Syringe rack (Ref.: 610300) or other syringe holder

(D) = Disposable
(R) = Reusable