



Annual Preventative Maintenance & Extended Service

As a Pro-Nox™ user, CAREstream America offers Annual Preventative Maintenance and Extended Service programs for your device mixer. Some programs consist of loaner devices and the more advanced programs include all Parts and Labor for necessary device repairs on top of an Annual Maintenance. The Pro-Nox™ manufacturer recommends Preventative Maintenance every 6 months.

Annual Preventative Maintenances include:

- Testing the individual features of the PRO-NOX
- Oxygen and Nitrous Oxide concentration
- Low input pressure alarms and shutoff functions
- 87 PSI shock and leak test
- Demand valve sensitivity and exhalation resistance
- Emergency air intake resistance
- Vibration
- Calibration & Adjustments needed to bring back into compliance

Preventative Maintenance Program	Price Per Year*	Loaner Included	Repairs Included
PNS-ANLPM – ANNUAL PREVENTATIVE MAINTENANCE	\$690.00	No	No
PNS-ANLPML – ANNUAL PREVENTATIVE MAINTENANCE w/loaner	\$1,190.00	Yes	No

Extended Service Contract benefits are:

- Fixed maintenance cost that allows accurate yearly budgeting
- ALL parts, labor and shipping INCLUDED
- Less demand on your technical personnel and resources
- Includes two annual preventive maintenance

Service Contract Program	Price Per Year*	Loaner Included	Repairs Included
PNS-EXTSCT – EXTENDED SERVICE CONTRACT	\$1,190.00	No	Yes
PNS-EXTSCTL – EXTENDED SERVICE CONTRACT w/loaner	\$1,690.00	Yes	Yes

**Taxes additional, where applicable*

