

STERILOCK™ MICROSYSTEM

Microneedling for the Modern Provider

Designed by skincare professionals for skincare professionals. The SkinStylus® SteriLock™ MicroSystem provides superior technology with multiple cartridge arrays at your fingertips for safe and exceptional results.

First introduced nearly a decade ago, SkinStylus has re-imagined microneedling with the 4th generation of their FDA approved SteriLock™ MicroSystem. Featuring both corded and battery operated models, our plug-and-play device offers interchangeable components to customize SteriLock™ to best suit the needs of you and your patients.

The SteriLock™ MicroSystem is FDA cleared up to 2.5 mm and is the only approved device for sheath-free operation. Other systems can trap biohazardous fluids inside the device that cannot be removed or reprocessed in an autoclave. The SteriLock™ is the only microneedling device with a removable battery and a nosecone that can either be autoclave sterilized or intermediate disinfected. Ask about our dependable, multi-year warranty and available trade-in program*.

The SkinStylus® Difference

- Corded & battery powered (long battery life)
- Patented Nosecone for sterile operation & device protection; NO sheath required
- FDA cleared Class II for atrophic scars
- Built-in redundancy & multi-year warranty
- Variety of cartridge arrays, including a novel 36-pin Hi-Lo option
- 3x more microchannels per second than closest competitor



Safety & Security with the Biolock® Cartridge

The patent-pending Biolock cartridge is made from 308 18/8 medical grade steel with FDA tested Makrolon polycarbonate. Our durable Biolock cartridges ensures patient safety through a 7-step system aimed to prevent leaks and contamination. Cartridges include a safety cap and are offered in 3 different configurations for tailored microchannels: 12 pin, 36 pin and 36 Hi-Lo pin.




Illustration of the interior actions of a Biolock cartridge.


855.892.3872 @ info@carestreamamerica.com carestreamamerica.com

*Trade-in conditions apply. Product distributed by CAREstream America



FDA APPROVED STERILOCK™ MICRONEEDLING SYSTEM

Our SkinStylus SteriLock® MicroSystem has been cleared by the FDA as a Class II Medical Device for the indication of treating abdominal scars (such as tummy tuck scars) to 2.5 mm, in adults over 22 years old. (510k #K200044) SkinStylus® MicroSystems provide advanced technology that offers exceptional results.



Vertical hyperpigmentation mark and horizontal stretch marks on 57 year old female. After 14 days, client was ecstatic to see the dark mark gone and an improvement in the appearance of the skin around the area.

A MODULAR DESIGN FOR THE MODERN PROVIDER

Two models available with accessories to ensure zero downtime between patients: Quickstart Kit or Professional Kit



| QUICKSTART KIT | PROFESSIONAL KIT |
|------------------------------|------------------------------|
| One Pen - Corded Only | Two Pens - Corded & Cordless |
| ABS Case + Acrylic Stand | ABS Case + Acrylic Stand |
| Motor | Motor |
| Nosecones (2) | Nosecones (4) |
| Corded Speed Head Controller | Corded Speed Head Controller |
| Modular Li-on Batteries (2) | Modular Li-on Batteries (2) |
| 9' Power Cord & AC Block | 9' Power Cord & AC Block |
| 3' Charging Cord & AC Block | 3' Charging Cord & AC Block |
| Cartridges (5) | Cartridges (10) |
| 2-Year Warranty | 3-Year Warranty |
| Patient Brochures (50) | Patient Brochures (50) |

CONNECT TO SCHEDULE A DEMO!

855.892.3872
 info@carestreamamerica.com
carestreamamerica.com

Product distributed by CAREstream America

