CAREROVERTM + ADAPTER ARM

HD vein imaging meets convenience with the **CARErover Rolling Stand** & Adapter Arm bundle

The near-infrared (NIR) technology of VeinViewer® Flex is now more portable than ever by adding the CARErover and Adapter Arm bundle. The CARErover utility cart was designed specifically to integrate with the VeinViewer® Flex system. The Adapter Arm flex stud easily connects to the existing VeinViewer® Flex Arm Mount to maintain motility when repositioning the vascular access device. Accuracy, stability and mobility combine for the ideal vein imaging package.

CARErover HEIGHT 43.5" (with Adapter Arm) Approx. 32" (with Adapter Arm & ADAPTER ARM REACH VeinViewer® Flex System) **INCLUDED TOOLS** 3/16" Allen Key & 1/2" Socket Wrench **ACCESSORIES** Transport handle & storage basket

Combine Pro-Nox and VeinViewer® Flex Systems by adding an Adapter Arm to your existing rolling stand.

VeinViewer® Flex VeinViewer System with Arm Mount (not included) Adapter Arm with flex stud to reposition VeinViewer® Flex

> Adjustable height handle for comfort & easy maneuverability

Storage basket for additional supplies & resources

39" steel stability pole

18" diameter, 5-caster star base with locking wheels for smooth mobility

Connect to learn more or schedule a private demo.









(carestreamamerica.com

